



# Technical data sheet



	Oticon Opn 1	Oticon Opn 2	Oticon Opn 3	
<b>Speech Understanding</b>	OpenSound Navigator™	Level 1	Level 2	Level 3
	Balancing power effect	100%	50%	50%
	Max. noise removal	9 dB	5 dB	3 dB
	Speech Guard™ LX	Level 1	Level 2	Level 3
	Spatial Sound™ LX	4 estimators	2 estimators	2 estimators
	Soft Speech Booster LX	•	•	•
	Speech Rescue™ LX	•	•	•
<b>Sound Quality</b>	Clear Dynamics	•	•	-
	Spatial Noise Management	•	•	-
	Fitting Bandwidth*	10 KHz	8 KHz	8 KHz
	Processing Channels	64	48	48
	Bass Boost (streaming)	•	•	•
<b>Listening Comfort</b>	Transient Noise Management	4 configurations	On/Off	On/Off
	Feedback shield LX	•	•	•
	Wind Noise Management	•	•	•
	Binaural Coordination***	•	•	•
<b>Personalization &amp; Optimizing Fitting</b>	YouMatic™ LX	3 configurations	2 configurations	1 configuration
	Fitting Bands	16	14	12
	Multiple Directionality Options	•	•	•
	Adaptation Management	•	•	•
	Oticon Firmware Updater	•	•	•
	Fitting Formulas	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0
	Acoustic Notifications	•	•	•
<b>Connecting to the World</b>	Stereo streaming (2.4 GHz)	○	○	○
	Oticon ON App	○	○	○
	ConnectClip	○	○	○
	Remote Control 3.0	○	○	○
	TV Adapter 3.0	○	○	○
	Autophone	○	○	○
	Tinnitus SoundSupport™***	•	•	•
<b>Battery life, hours**</b>	55-60 / 105-115	55-60 / 105-115	55-60 / 105-115	

\* Bandwidth accessible for gain adjustments during fitting  
 \*\* Battery size 312 - IEC PR41 / Battery size 13 - IEC PR48.  
 Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).  
 \*\*\* If push button is chosen  
 • Default ○ Optional - Not included

## OTICON | Opn ITC, ITE HS & FS 90



Oticon Opn™ ITC, ITE HS & FS introduce an updated faceplate design.


OpenSound Navigator™ provides better speech understanding by continuously analyzing the environment, balancing all sound sources and attenuating the dominating noise.

TwinLink™ wireless technology combines binaural communication and 2.4 GHz connectivity in stereo directly to external digital devices with very low power consumption. 2.4 GHz is an optional.

Oticon Opn is a Made for iPhone® hearing aid.

Oticon Opn is built on the Velox™ platform, providing frequency resolution in 64 channels (Opn 1).

Fully programmable with updatable firmware, the Velox platform is ready for the future.



For information on compatibility, please visit [www.oticon.com](http://www.oticon.com)

Call  
 Absolute Hearing Solutions  
 For Hearing Aid Deals  
 614-452-4280  
[www.absolutehearingsolutions.com](http://www.absolutehearingsolutions.com)




Technical data Measured according to		Ear Simulator IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010			ZCC Coupler ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006		
Oticon Opn ITC ITE HS & FS 90		Opn 1	Opn 2	Opn 3	Opn 1	Opn 2	Opn 3
Frequency range Hz		110-9500	110-7500	110-7500	100-7900	100-7500	100-7500
OSPL90	Peak	131 dB SPL			120 dB SPL		
	1600 Hz	123 dB SPL			115 dB SPL		
	HFA-OSPL90	124 dB SPL			116 dB SPL		
Full-on gain*	Peak	65 dB			55 dB		
	1600 Hz	57 dB			48 dB		
	HFA-FOG	58 dB			50 dB		
Reference test gain		48 dB			39 dB		
Telecoil output (1600 Hz)	1 mA/m field	87 dB SPL			-		
	10 mA/m field	107 dB SPL			-		
	SPLITS L/R	-			96/96 dB SPL		
Total harmonic distortion (Input 70 dB SPL)	500 Hz	2 %			< 2 %		
	800 Hz	2 %			< 2 %		
	1600 Hz	2 %			< 2 %		
Equivalent input noise level	Omni	18 dB SPL			15 dB SPL		
	Dir	28 dB SPL			27 dB SPL		
Battery consumption**	Typical	1.8 mA			1.8 mA		
	Quiescent	1.7 mA			1.7 mA		
Battery life, calculated, hours 312 and 13***		100 / 175			100 / 170		
IRIL (IEC 60118-13:2016)		700/1400/2000 MHz: 20/12/6 dB SPL					

\* Measured with the gain control of the hearing aid set to its full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.

\*\* Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.

\*\*\* Based on the standardised battery consumption measurement (IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.



<b>Operating conditions</b> Temperature: +1°C to +40°C  Relative humidity: 5% to 93%, non-condensing	<b>Storage and transportation conditions</b> Temperature and humidity should not exceed the following limits for extended periods during transportation and storage.  Temperature: -25°C to +60°C Relative humidity: 5% to 93%, non-condensing
------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.

Call  
 Absolute Hearing Solutions  
 For Hearing Aid Deals  
 614-452-4280  
[www.absolutehearingsolutions.com](http://www.absolutehearingsolutions.com)

**oticon**  
 PEOPLE FIRST