

ReSound LiNX 3D Rechargeable Battery Care Tips



Deep discharging of the rechargeable battery can significantly reduce the capacity of the battery to hold a full charge. The rechargeable batteries can recover from this but it can take 5 to 10 charging cycles. You can prolong the life of your rechargeable battery with the following guidance:

1. Do not reboot hearing instruments after low battery warning beeps.

When the low battery warning is heard, it is an indication to immediately charge the hearing instruments in the charger or replace with zinc air batteries.

2. Fully charge the hearing instruments every day.

It can take up to 7 hours for the batteries to be fully charged depending on how much charge is left in the batteries.

3. Verify that the hearing instruments are charging.

The green light on the charger will start flashing immediately when hearing instruments are placed in the charger indicating the charging process has started.

If a solid red light appears, the hearing aids are not charging. Ensure there are rechargeable batteries in the hearing instruments and that they are properly inserted into the charger.

4. Keep the hearing instruments in the charger when not in use.

Do not leave the hearing instruments with the battery doors open

5. Remove the batteries when using drying devices.

The temperature ranges of some drying devices may exceed the battery's operating temperature damaging the batteries.



Absolute Hearing Solutions LLC
absolutehearing@att.net

750 Cross Pointe Road Suite F
Gahanna, Ohio 43230

Phone: 614-452-4280
Toll Free: 888-803-2159
Fax: 614-577-0481

www.absolutehearingsolutions.com

Manufacturer according to FDA:

GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-888-735-4327
resound.com

Manufacturer according to Health Canada:

ReSound Canada
303 Supertest Road
Toronto, Ontario M3J 2M4
Canada
1-888-737-6863
resound.com



MK605089 Rev A 2018.05

ReSound GN