

propel™

WIRELESS HEARING PRODUCTS



Call
Absolute Hearing Solutions
For Hearing Aid Deals
614-452-4280
www.absolutehearingsolutions.com

micro tech



the carefree way to
**CONNECT TO
YOUR WORLD**

Propel hearing aids by MicroTech are designed to make listening easy and enjoyable again. With their advanced features, they're engineered to:

- Reduce listening effort in noisy environments*
- Improve listening to music and TV
- Enhance phone calls
- Virtually eliminate whistling and buzzing
- Fit comfortably, customized for your ears



Propel hearing aids feature MicroTech's latest noise reduction and speech preservation system, which is designed to deliver incredible sound clarity even in the noisiest situations.

Propel can also stream stereo sound from your TV, radio, computer and even your cell phone directly to your hearing aids. With Propel, you can enjoy hands-free phone conversations and listen to TV or radio at the volume you want, while others around you can listen at the volume that's most comfortable for them.

If you have been looking for a hearing aid designed to make both conversations comfortable, and TV, music and phone conversations sound like never before, Propel is it.

Life

without limitations



Combined with our innovative SurfLink® accessories, Propel hearing aids can also make it easier to enjoy TV, music, cell phone conversations and more — no wires attached!

SurfLink Media

This set-and-forget media streaming solution connects to TVs, MP3 players and more to wirelessly stream audio directly to your Propel hearing aids, so you can hear the volume at the level you want, while everyone else in the room can listen at the volume they want.



SurfLink Remote

This handy remote lets you adjust and control your Propel hearing aids with the push of a button.



SurfLink Mobile

Enjoy hands-free phone conversations with this one-of-a-kind device. That's because SurfLink Mobile turns your Propel hearing aids into both your cell phone microphone *and* receiver:



Microphone — Your hearing aids can pick up your voice and send it to the person you're talking to

Receiver — SurfLink Mobile streams the voice on the other end of the phone directly to both of your hearing aids

But that's just the beginning. SurfLink Mobile is not only a cell phone transmitter, it's also an assistive listening device, media streamer, and hearing aid remote all rolled into one.

propel™
w50
Vibrant

designed to
fit your life.

LIFESTYLES

CROWD



CAR



THEATER



PARTY



INDUSTRIAL
SOUNDS



MALL



MEETING



RESTAURANT



OUTDOORS



WORSHIP



ONE-ON-ONE



HOME



propel™
w40

Active

No matter where you are or what you're doing, Propel is designed to make listening easier. It's the ideal hearing aid for people who are active and on the go, enjoy spending relaxing time at home — or both.

propel™
w30

Social

Binaural Spatial Mapping

PREMIUM

ADVANCED

WIRELESS
FEATURES

Media Connectivity



Synchronized User Controls



Cell Phone Connectivity





engineered to make
HEARING CAREFREE

Technology

Benefits

Strategic Feedback Control

Virtually eliminates buzzing and whistling

SpeechTrax²

Designed to maintain speech understanding in noisy environments and reduce listening effort*

Dynamic Directional Detection

Enhances SpeechTrax² by pinpointing sounds you want to hear and minimizing those you don't

Phone

Used with SurfLink Mobile, can enable hands-free cell phone conversation, while Automatic Telephone Response detects when you're on the phone and adjusts setting for optimal listening

Music & Television

Automatically adjusts to the preset music and TV listening settings most comfortable for you — while optional SurfLink Media can stream audio from TV, computer or MP3 player directly to your hearing aids

HydraShield^{®2}

Makes your hearing aids resistant to water, wax, sweat, oil and corrosion

FrequencyTrax

Replicates hard-to-hear high frequency sounds in lower frequencies, designed to enhance and expand audibility

Self Learning

Gradually and automatically learns your volume control preferences for less manual adjustments

Convenient Intuitive Features

Automatic voice and tone indicators tell you the status of your hearing aid and your hearing aid battery

Ear-to-Ear Wireless Communication

Coordinates communication between hearing aids, eliminating the need to manually adjust both devices

Features may vary by hearing aid style or technology level



propel™

STYLES and COLORS


Easy-to-adjust Propel is available in two styles and a variety of discreet and stylish colors. Consult your hearing care professional and start experiencing Propel for yourself.

mini Behind-The-Ear (mini BTE)




Available July 2012

Standard Colors



Bright Colors



Receiver-In-Canal (RIC)





Absolute Hearing Solutions

Absolute Hearing Solutions LLC

absolutehearing@att.net

750 Cross Pointe Road Suite F
Gahanna, Ohio 43230

Phone: 614-452-4280
Toll Free: 888-803-2159
Fax: 614-577-0481

www.absolutehearingsolutions.com



6425 Flying Cloud Drive
Eden Prairie, MN 55344

800.745.4327
www.microtechhearing.com

Individual results may vary.

*Sarampalis, A., Kalluri, S., Edwards, B., Hafer, E. (2009, October). Objective measures of listening effort: Effects of background noise and noise reduction. *Journal of Speech, Language, and Hearing Research*, 52, 1230-1240.