





























HEARING AID STYLE COMPARISON CHART


 STYLE	 FULL NAME	 VISIBILITY	 RECHARGEABLE	 BLUETOOTH	 BEST FOR
IIC 	Invisible-In-Canal	 Invisible	 No	Limited	Cosmetic appeal
CIC 	Completely-In-Canal	 Very discreet	 Some models	 Some models	Mild hearing loss
ITC 	In-The-Canal	 Small	 Yes	 Yes	Easy handling
ITE 	In-The-Ear	 Moderate visibility	 Yes	 Yes	Dexterity & controls
RIC 	Receiver-In-Canal	 Discreet behind-the-ear	 Yes	 Yes	Natural sound & mild to severe hearing loss
BTE 	Behind-The-Ear	 Most visible	 Yes	 Yes	Mild to profound hearing loss



**BETTER HEARING STARTS HERE.
BETTER CALL TODAY!**

Contact your Hearing Experts for more details at

ABSOLUTE HEARING SOLUTIONS

 750 Cross Pointe Road, Suite F
Columbus, Ohio 43230

 614-452-4280

 www.absolutehearingsolutions.com