# G E N E S I S<sup>AI</sup>

# Wireless Accessories & Mobile Apps

Call Absolute Hearing Solutions For Hearing Aid Deals 614-452-4280 www.absolutehearingsolutions.com





# Stay connected to the world.

Our versatile line of **2.4 GHz StarLink accessories** connects wireless Genesis AI hearing aids to most modern smartphones, televisions or media devices so patients can easily stream phone calls, music, voice assistants and more.





### **Table Microphone**

- Enhances listening in group settings like meetings, family gatherings or noisy restaurant environments.
- Place it in the center of a table and the eight built-in microphones will home in to the primary speaker's voice and stream it directly to the Genesis AI hearing aids.

### **TV Streamer**

- Streams audio directly from TVs or other electronic audio source directly to Genesis Al hearing aids.
- Easy to use and supports both analog and digital input sources.



### **Remote Control**

- Lets user control memory and volume, mute hearing aids, and turn other special features on and off.
- Ideal for patients with older mobile devices.



### Remote Microphone +

- Provides extra listening assistance in noisy environments to hear one-on-one conversations.
- Easily streams binaurally from a variety of audio sources, including Bluetooth®, Telecoil, FM or Line-In.



### **Mini Remote Microphone**

- Enhance one-on-one conversations in noisy environments with this small remote mic.
- Clips onto the clothing of the person being spoken to.
- Can be used as a TV streamer by placing near the sound source.



Download the My Starkey or Hear Share app here

### ALL-NEW My Starkey app

Designed with the input of hearing aid wearers, My Starkey easily pairs with Genesis Al hearing aids to help patients enjoy greater listening moments, improved wellbeing, and easier daily living. All from the palm of their hand. It's optimized for readability and places the most-used features right where patients need them.



### My Starkey user-friendly features:

ß

### Hearing control

Adjust hearing aid settings, stream calls and more.

#### Personal health monitoring

Count steps, track engagement and achieve health goals.

#### 🚡 Translate

Translate 77 languages in real time directly through the hearing aids.

### Edge Mode+ Uses Al for

oses Ai for on-demand sound quality optimization.

## Binaural streaming

Enjoy real, not fake, binaural streaming.

#### Smart Assistant

S

Voice activated audio adjustments, reminders, answers to hearing aid topics, and more.

### Fall detection

Detect falls and alert family or friends.

#### Patient education


A complete library of fresh "How-to" videos.



Scheduled tasks, appointments, and more can be streamed right to the hearing aids.

9:41	¢lis.	-	
ane			
🕒 Left 🗈 Rig	9:41		- + + h.
< Jan 28	3, 202		
B thr Daily goa	il is 12	6	
Steps 7,000 Daily goa	l is 10,0	<b>EAR</b> Shar	e
Standing Standing Shr Daily goa			
(in the second s	) mir		

### Thrive Care is now Hear Share



wearers share select progress markers with friends and family members, including:

- Their latest strides in hearing aid use.
- Physical activity and social engagement goals.
- Fall alerts and notifications.





#### **Absolute Hearing Solutions LLC** absolutehearing@att.net

750 Cross Pointe Road Suite F Gahanna, Ohio 43230

Phone: 614-452-4280 Toll Free: 888-803-2159 Fax: 614-577-0481

www.absolutehearingsolutions.com

Starkey logo, Genesis, Edge Mode, Hear Share and StarLink are trademarks of Starkey Laboratories, Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Starkey is under license. ©2023 Starkey Laboratories, Inc. All Rights Reserved. 1/23 BKLT3176-00-EN-ST